

PORTER

SENTINEL SYSTEM™ MANIFOLD

DO NOT CHANGE PIPE SIZE

DO NOT CHANGE PIPE SIZE

NITROUS OXIDE
PIPELINE 3/8"
O.D. COPPER

OXYGEN PIPELINE
1/2" O.D. COPPER

- 1) USE TYPE K OR L, PRE-CLEANED, DEGREASED, CAPPED COPPER ONLY.
- 2) USE SILVER SOLDER OR SIMILAR BRAZING ALLOY WITH AT LEAST A 1000°F MELTING POINT ONLY.
- 3) TEST FOR LEAKS WITH DRY NITROGEN AT 150 PSI FOR 24 HOURS.

ACTUAL SIZE TEMPLATE

