



NOTES:

1. MATERIAL - 5MIL CLEAR MYLAR.
2. BLACK GRAPHICS AND LETTERING

-	ECO-0187634	ECO RELEASE FOR PRODUCTION	A. MARQUES	FEB. 14, 2017	
REV	EC NO.	DESCRIPTION	NAME	DATE	APPVD
REVISIONS					
				FILENAME: LB-730	
DATE	DR. A. MARQUES	CH.	APP.	SCALE 1:2	NUMBER
OCT. 10, 2016					LB-730
					SH 1 OF 1
					REV
					-

**PROPRIETARY NOTICE**  
THIS DOCUMENT AND THE DATA DISCLOSED HEREIN IS NOT TO BE REPRODUCED, USED OR DISCLOSED WHOLE OR IN PART TO ANYONE WITHOUT THE PERMISSION OF PORTER INSTRUMENT DIVISION, PARKER HANNIFIN CORPORATION.



Parker Hannifin Corporation  
Porter Instrument Division  
245 Township Line Road  
Hatfield, PA 19440, USA

TEMPLATE, HOLE LOCATION,  
SENTRY SEDATE SCAV

LB-730  
SH 1 OF 1

LB-730  
REV -

MOUNTING HOLE LOCATION

ALIGN WITH COVER DETAIL

OUTER EDGE REFERENCE

LB-730  
REV -