

Porter Nitronox™ Field Unit

Portable, Demand Flow Nitrous Oxide
And Oxygen Analgesia System

The Trusted Name In Nitrous Oxide

Patient Self-Administered
Fixed 50/50 N₂O and O₂
Fast Acting
Short Duration Of Effect



The Porter Nitronox™ Field Unit can provide on-the-spot inhalation analgesia, bringing relief to suffering patients at the scene of an injury.



Nitronox is a patient self-administered, demand-flow system that delivers a safe and effective mixture of 50% nitrous oxide and 50% oxygen through a hand-held demand valve. The Nitronox system easily connects to small, custom nitrous oxide cylinders and your standard oxygen supply.

Pain relief starts within minutes after inhalation. Once inhalations are discontinued, the effect of the nitrous oxide is quickly reversed by breathing room air.

That means you can take care of your patients, while they manage their own pain safely and effectively.

Nitronox – For effective and fast, yet safe pain relief.

Get more information at: www.porterinstrument.com/nitronox

The Porter Nitronox™ Field Unit provides fast, effective and convenient pain management in emergency settings.

Control your patient's pain and anxiety at the scene of an injury, with the fast, effective Porter Nitronox™ Field Unit, a compact nitrous oxide system that fits easily into an ambulance.

The features of the Nitronox unit protect your patients while relieving pain.

- **Patient self-administered**

A hand-held supply valve is activated only by the patient's inhalation.

- **Preset 50/50 mixture of N₂O/O₂** eliminates the risk of delivering a hypoxic dose.

- **100% O₂ capability** In the event of N₂O depletion, the patient gets 100% oxygen.

- **O₂ enrichment** In the event of shallow breathing, the unit will automatically increase O₂ percentage.

- **Quick disconnect** Allows for easy removal and storage of the patient demand valve.



NOX-5000E Nitronox Field Unit Package

Includes Nitronox Mixer, Nitrous Oxide Demand Valve (with delivery hose and quick connect), Oxygen Hose (with DISS connection), 2 Nitrous Oxide Cylinders (empty), and Soft Carry Case.

Service and Calibration Testing available.

Call for information on Trade-In Options.



Contact Information:

Porter Instrument Division
Mike Civitello – Medical Product Manager
Email: michael.civitello@parker.com
Tel: 215-660-8224
Fax: 215-723-5106

www.porternitronox.com

Exclusive U.S. Distributor:

 **HENRY SCHEIN®**
MEDICAL | EMS

Henry Schein, EMS
Tel: 1-800-845-3550
Fax: 1-800-533-4793
Email: scott.bruner@henryschein.com



• **SAFE** • **PATIENT SELF-ADMINISTERED**
• **FAST-ACTING** • **SHORT DURATION OF EFFECT**