


# Porter **Sentry**Sedate<sup>HP</sup>

## Portable 4 Cylinder Nitrous Oxide and Oxygen System



*Porter MXR-D Digital and Analog Flowmeter*



*Porter MXR-1 Analog Flowmeter*

Porter's newest Nitrous Oxide Sedation System features portability and security with:

- Dual Locking Doors
- Holds 2 Oxygen and 2 Nitrous Oxide "E" Cylinders
- Oxygen/Nitrous Oxide Quick Connect Ports
- Internal and External Storage
- Working Surface
- Durable Locking Wheels

### **Models**

SENTRY-HD34  
MXR-1 Hospital Sentry Sedate  
SENTRY-HD44  
MXR-D Hospital Sentry Sedate

### **Dimensions & Weight**

Height - 33"  
Depth - 24"  
Width - 18"  
Weight - 96 lbs (add 2 lbs for Porter Digital MXR) - does not include tanks



[www.porterinstrument.com/medical](http://www.porterinstrument.com/medical)

# Sentry Sedate<sup>HP</sup> Features For Portability, Security And Convenience



The sturdy cabinet design features both internal and external storage.



Sentry Sedate holds two nitrous oxide cylinders and two oxygen cylinders.



A side hook conveniently holds the breathing circuit and hoses.



Large locking wheels allow the unit to be easily moved from room to room.



The dual locks ensure safe, secure cylinder storage.



Oxygen and nitrous oxide quick connect ports bypass cylinders when wall gas is available.

**For more information, contact Praxair - Authorized Porter Distributor:**

Ph: 1-877-798-0854 Email: [infohealth@praxair.com](mailto:infohealth@praxair.com)

[www.PorterInstrument.com/medical](http://www.PorterInstrument.com/medical)

