

# DIAGRAM, YOKE CHECK VALVE REPAIR KIT

---

---

# DIAGRAM, YOKE CHECK VALVE REPAIR KIT

QTY.		PARTS LIST		
-001	-000	DESCRIPTION	PART NO.	ITEM
-	1	CHECK VLV.POPPET,O2,VT ELAST.	PA-616-000	1
1	-	CHECK VLV.POPPET,N2O,NITR ELAS	PA-616-001	
-	1	O2 YOKE CHECK VLV SEAT SCREW	PA-617-000	2
1	-	N2O YOKE CHECK VLV SEAT SCREW	PA-617-001	
1	1	3/32 HEX WRENCH	PA-623-1	3
1	1	5/32 HEX WRENCH	PA-623-3	4
1	1	SCREW,FILTER RETAINING	PA-768-000	5
1	1	DIAGRAM,YOKE CHK VLV REPAIR KT	FM-541	6

