



# PORTER **MIDAS**<sup>TM</sup>

*Digital Touchscreen Flowmeter*



**Innovation and  
technology  
at your fingertip**



# PORTER

*The Trusted Name In Nitrous*

[www.PorterInstrument.com/Medical](http://www.PorterInstrument.com/Medical)

# Introducing Porter MIDAS™ – A new nitrous oxide flowmeter with all the right touches.



## Features & Benefits

With its advanced touchscreen interface, Porter Midas™ redefines precision and accuracy in nitrous oxide and oxygen delivery. Touch and Slide™ gas flow and concentration control *and so much more with a touch of the finger!*

- Touch and Slide™ flow and concentration controls\*
- High flow capability – up to 18 LPM
- Capacitive touch
- Ergonomic and intuitive controls
- Scavenger Plus Interface – Effective scavenging when using disposable breathing circuits
- Optional eAVS – Electronic Automatic Vacuum Switch – precision scavenger control\*
- Audible and visual alerts
- Porter Midas™ App and iPad™ exportable procedure log
- Porter Midas™ App and iPad™ remote control option
- Gas used and gas remaining log
- End-of-procedure O<sub>2</sub> recovery program
- 2-year warranty

## Safety Features

- Oxygen fail-safe
- Optional nitrous lock-out PIN for security
- Emergency air intake valve
- Diameter Indexed Safety System connections
- Oxygen flush
- Non-rebreathing check valve
- Auxiliary oxygen port for demand valve or oxygen flowmeter
- Self-check diagnostic test during start-up and ongoing during use

\* eAVS and Touch and Slide™ are patent pending

# Scavenger Plus




- Allows connection to wall source VAC / WAGD outlets for scavenging of exhaled waste gas
- Visual scavenger control gauge for precise control of vacuum flow
- Exhalation bag for visual confirmation of adequate scavenging flow
- Quiet operation
- Required if using full facemask or mouthpiece

## eAVS: Electronic Automatic Vacuum Switch



- Digital scavenging display and control on Porter Midas™ touchscreen
- Precise control of vacuum flow
- Automatically activates vacuum flow when gas flow begins
- Required if using Silhouette nasal hoods



Visit our website for more details: [www.PorterInstrument.com/Medical](http://www.PorterInstrument.com/Medical)

# Porter MIDAS™ Mounting Options



**2-Cylinder Mobile Cart**




**Wall Mount**

## Delivery Options



**Disposable Breathing  
Circuit and Mask**



**Silhouette Nasal Hood**

**PORTER**  
**MIDAS™**  
*Digital Touchscreen Flowmeter*

**Contact Porter for more information:**

P: 215-723-4000 Email: [porternitrous@parker.com](mailto:porternitrous@parker.com)

[www.PorterInstrument.com/Medical](http://www.PorterInstrument.com/Medical)

**PORTER**  
*The Trusted Name In Nitrous*

