

PORTER

The Trusted Name In Nitrous

Nitrous Oxide and Oxygen for Pain and Anxiety Management During In-Office Gynecologic Procedures



Work More Efficiently and Increase Treatment Acceptance

Nitrous oxide and oxygen ($\text{N}_2\text{O}/\text{O}_2$) provides a safe, simple and fast-acting alternative to oral medications for minimal sedation. During the procedure, patients experience reduction in pain and anxiety due to the analgesic and anxiolytic properties of N_2O . After the procedure is completed, the effects on the patient wear off within minutes. OB/Gyn procedures where $\text{N}_2\text{O}/\text{O}_2$ is utilized include:

- IUD insertion
- Pelvic exams
- Endometrial ablation
- Hysteroscopy
- Labiaplasty
- Biopsies
- LEEP procedure
- Vaginal rejuvenation

NITRONOX PLUS®

Demand Flow Nitrous Oxide & Oxygen Systems

The next generation of Nitronox® features a sleek new design and a choice of fixed or adjustable nitrous oxide, PLUS so much more! The Nitronox® name is synonymous with quality and the ability to effectively reduce patient pain and anxiety during minimally invasive medical procedures. Nitronox Plus® builds on that success.

Porter has taken the benefits of an adjustable continuous flow system and combined them into a new demand flow system that allows for a predictable patient experience. Now available are options that offer 0-50% or 0-70% adjustable nitrous oxide along with a new 50/50 fixed system. These 3 new models allow you to choose the system that best fits your needs. Regardless of application, Nitronox Plus has the right option for you.



50/50 Fixed O₂ / N₂O



0-50% Adjustable N₂O



0-70% Adjustable N₂O



Nitronox Plus on 2-Cylinder Stand

Nitronox Plus® – Demand Flow with enhanced features

- Available in fixed or adjustable mixture concentration models
- Key lock mechanism for safety and security
- Audible low gas pressure alarm
- Internal, but accessible demand valve
- Safety function test button

Designed to meet infection control standards

- One-way directional flow of inhaled and exhaled gas
- Single-use, disposable breathing circuits
- Single-use masks and mouthpieces
- Front panel user interface prevents surface cleaners and disinfectants from entering the device



PORTER

The Trusted Name In Nitrous

P: 888-723-4001 E: porternitrous@parker.com



www.PorterInstrument.com/obgyn

FM #1523