

Porter Nitronox Plus Nitrous Oxide Analgesia System: Empower Your Patients With Nitrous Oxide!



The Nitronox® name is synonymous with quality and the ability to effectively reduce patient pain and anxiety during labor and delivery. Nitronox Plus builds on decades of that success and is the culmination of years of gathering user feedback to design a device that offers healthcare professionals the functionality they've asked for.

Equipment

- Nitronox Plus is a demand flow system with fixed 50/50 nitrous oxide and oxygen
- Key lock system for safety and security
- Low gas pressure alarms alert you when to change cylinders
- Internal, but accessible demand valve limits downtime during preventive maintenance checks
- Visual gauges to confirm device is ready for use
- Available with multiple cart and wall-mount options



Patients

- Manage pain and anxiety during labor and post-birth repairs
- Patient self-administered with rapid onset of effects
- Attract new patients and improve patient satisfaction

Porter Instrument

- Leading manufacturer of nitrous oxide and oxygen delivery systems for over 50 years
- Resources include: Sample policies, contacts at other hospitals, in-service training, online training and more
- Manufactured and serviced in the USA

PORTER
The Trusted Name In Nitrous



SAFE • PATIENT SELF-ADMINISTERED
• FAST-ACTING • SHORT DURATION OF EFFECT

Nitronox Plus Nitrous Oxide and Oxygen Analgesia System

What Patients Are Saying

Here are some quotes from a recent survey showing what new mothers had to say about using Nitronox during and after labor.

"Knowing I would have access to a pain aide that would not interfere with the health of my baby or my ability to move and labor in the way I wanted was very appealing."

"I liked that I was in total control and was able to regulate how much or how little I needed. I also liked the fact that the gas did not take away my ability to feel the pain or listen to what my body was telling me to do."

"During stitching after the birth - it caused me to relax and not feel the anxiety I usually have with needles or any pain associated with the experience."

Complete Scavenging Solution with Scavenger Plus

- Connects to all wall vacuum and WAGD outlets to scavenge exhaled waste gas
- Adjustable vacuum control with gauge to ensure proper setting
- Breathing bag offers visual confirmation of gas being scavenged
- Quiet operation



Visit our website to see articles and videos about the benefits of Porter Nitronox Plus!
www.porterinstrument.com/labor

PORTER
The Trusted Name In Nitrous

www.PorterInstrument.com/Labor

Parker | Porter

SAFE • PATIENT SELF-ADMINISTERED
• FAST-ACTING • SHORT DURATION OF EFFECT

Email: porternitrous@parker.com